

# photomultiplier amplifier-discriminator

## AD8 data sheet

### 1 description

The AD8 amplifier-discriminator is a compact, low noise, high gain electronics module designed to be used in photomultiplier photon counting systems to generate TTL output pulses. An output pulse of fixed width and amplitude is produced for each photomultiplier output pulse that exceeds the input discriminator threshold level. The threshold level is fixed to -2mV, which is sufficiently sensitive to enable most photomultipliers with 10 or more dynodes to operate in single photon counting mode.

Two outputs are provided; one with a pulse-pair resolution of 25ns for input count rates up to 150Mcps, and the other with a pulse-pair resolution of 10 $\mu$ s for rates up to 100kcps. The slower output is intended to be used where it is required to eliminate photomultiplier afterpulses.

The AD8 is housed in a metal box and is compatible with the ET Enterprises MCS-CT3 multi-channel scaler/counter-timer. When used together, the AD8 can be powered from the MCS-CT3. These units are also compatible with the ET Enterprises range of ambient temperature and cooled photomultiplier housings, enabling high performance and cost-effective photon detection systems to be set up with minimal effort.

### 2 applications

- counting photomultiplier pulses exceeding a fixed threshold.
- single photon counting or applications involving multiple photodetection pulses.

### 3 features

- high input sensitivity and low noise
- count rates up to 150Mcps with dead time correction.
- signal and power compatible with the MCS-CT3 multi-channel scaler/counter-timer.
- input overload protection.
- can be used with negative or positive HV photomultiplier operation.



### 4 characteristics

	unit	min	typ	max
<b>input pulse amplitude</b>	mV	-1000		0
<b>discriminator level</b>	mV		-2	
<b>input impedance</b>	$\Omega$		50	
<b>input protection:</b> <b>limiting diodes</b>				
<b>supply voltage</b>	V	4.5	5	5.5
<b>supply current</b> <b>(@ 5V, no signal)</b>	mA		40	
<b>supply current</b> <b>(@ 5V, signal = 150Mcps)</b>	mA		110	
<b>output impedance</b>	$\Omega$		50	
<b>output pulse</b> TTL high level (terminated)	V	3.3		
<b>pulse-pair resolution:</b>				
17ns output	ns		25	
10 $\mu$ s output	$\mu$ s		10	
<b>dead-time</b>				
17ns output	ns		25	
10 $\mu$ s output	$\mu$ s		10	
<b>output pulse rise and fall time:</b>				
17ns output	ns		1.2	
10 $\mu$ s output	ns		1.2	
<b>signal rate count</b> <b>(without deadtime correction)</b>				
17ns output	cps			30M
10 $\mu$ s output	cps			100k
<b>signal count rate</b> <b>(with deadtime correction)</b>				
17ns output	cps			150M
<b>warm-up time</b>	s		1	
<b>temperature (operating)</b>	$^{\circ}$ C	5		55
<b>temperature (storage)</b>	$^{\circ}$ C	-40		60
<b>humidity (non-condensing)</b>	%			93
<b>weight</b>	g		160	


## 5 dynamic range

When using the fast (17ns) output, maximum dynamic range can be achieved by applying count rate correction to compensate for the departure from linearity due to pulse pile-up at high count rates. Dead time may be corrected for, as follows:

$$N = n/(1-nT)$$

where: N is the true count rate (cps),  
n is the measured count rate (cps),  
T is the count rate correction factor ( $25 \times 10^{-9}$ s),

Using this correction, deviation from linearity is typically within +/-5% at 120Mcps.




## 6 connections

The TTL output pulses from either the fast (17ns) or slow (10µs) BNC connectors may be connected to pulse counting equipment such as the ET Enterprises MCS-CT3 multi-channel scaler/counter-timer, or other compatible instrumentation. The AD8 can also be powered from the MCS-CT3, or from the auxiliary low voltage output of the ET Enterprises photomultiplier power supply HVLAB3000.

Although the AD8 is insensitive to most external sources of interference, it is recommended that the connection to the photomultiplier be made with a 50Ω terminated BNC screened cable of maximum length 50cm.

signal input	50Ω BNC socket
power input	2.1mm dc jack
signal output (17ns)	50Ω BNC socket
signal output (10µs)	50Ω BNC socket

## 7 external dimensions mm



## 8 ordering information

AD8	boxed electronic module
LEADAD8MCS-CT3	50 cm long power cable when used with MCS-CT3
AD8PSU	universal ac power adaptor
LEADBNC - 50cm	50 cm long 50Ω signal cable
LEADBNC - 100cm	100 cm long 50Ω signal cable

ET Enterprises Limited  
45 Riverside Way  
Uxbridge UB8 2YF  
United Kingdom  
tel: +44 (0) 1895 200880  
fax: +44 (0) 1895 270873  
e-mail: sales@et-enterprises.com  
web site: www.et-enterprises.com

ADIT Electron Tubes  
300 Crane Street  
Sweetwater TX 79556 USA  
tel: (325) 235 1418  
toll free: (800) 399 4557  
fax: (325) 235 2872  
e-mail: sales@electron tubes.com  
web site: www.electrontubes.com

an ISO 9001 and ISO 14001 registered company

The company reserves the right to modify these designs and specifications without notice. Developmental devices are intended for evaluation and no obligation is assumed for future manufacture. While every effort is made to ensure accuracy of published information the company cannot be held responsible for errors or consequences arising therefrom.

ET Enterprises  
electron tubes

© ET Enterprises Ltd 2015  
DS\_AD8 Issue 2 (30/03/15)